



**Product:** **Recombinant Human BMP-2**  
Powder

**Catalog #: 300-601P**  
**Lot #: G39L40A**  
**Assayed: Oct 2010**

**Date Received:** \_\_\_\_\_

**Expiration Date:\*** \_\_\_\_\_

| TEST                                 | METHODOLOGY   | RESULTS                    |
|--------------------------------------|---|----------------------------|
| Biological Testing<br>Endotoxin      | LAL   | 0.036 EU/μg                |
| Biochemical Testing<br>Purity        | SDS PAGE  | >95%                       |
| Biological Activity ED <sub>50</sub> | Alkaline Phosphatase activity induced in ATDC-5 Cells | 1.4 x 10 <sup>4</sup> U/mg |

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.*

*[Signature]* 19 Aug 2019

Created by/ date

*[Signature]* 19 Aug 2019

Verified by/ date

\*Expiration Date is one year from the date of receipt when stored at -20 to -80°C.