



**Product:** **Bovine Serum Albumin (BSA)**  
**Fraction V—Low Electrolyte**  
Powder

**Catalog #: 700-109P**  
**Lot #: C00J09K**  
**Origin: United States (100%)**

**Manufacture Date: Feb 2019**  
**Expiration Date: Feb 2022**

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Microbiological Testing			
Microbial Growth		Not Detected	Not Detected
Viral Testing	Modified 9 CFR 113.53c		
Extraneous Viruses (FAb)			
Bluetongue Virus		Not Detected	Not Detected
Vesicular Stomatitis Virus		Not Detected	Not Detected
Physical Testing			
pH @ RT, 10% Solution		4.9-5.2	5.1
Conductivity, 10% Solution		≤2,200 µS	370 µS
Biochemical Testing			
Purity	HPLC	≥98%	100%

*Filtered through a 0.2 µm filter before lyophilizing*

*This material was prepared from bovine plasma obtained from cattle of less than 30 months of age and that have received ante- and post-mortem health inspections from APHIS Inspectors within facilities inspected by the USDA.*

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.*

*Endy...*  
\_\_\_\_\_  
Created by/ date 09Apr2019

*[Signature]*  
\_\_\_\_\_  
Verified by/ date 09Apr2019