

Product: **Bovine Serum Albumin (BSA)**
Fraction V—Low IgG
 Powder

Catalog #: 700-105P

Lot #: C07H09H

Origin: United States (100%)


Manufacture Date: May 2016

Expiration Date: June 2019

<i>TEST</i>	<i>METHODOLOGY</i>	<i>SPECIFICATION</i>	<i>RESULTS</i>
Microbiological Testing			
Microbial Growth		≤3 cfu/mL	No Growth
Viral Testing			
Extraneous Viruses (FAb)			
Bluetongue Virus		Negative	Negative
Vesicular Stomatitis Virus		Negative	Negative
Physical Testing			
pH @ RT, 10% Solution		6.7-7.3	6.8
Conductivity, 10% Solution		≤3,000 μS	1050 μS
Biochemical Testing			
Purity	HPLC	≥98%	100%
Protease		Not detected	Not detected
IgG		Below detectable limits	Below detectable limits

Filtered through a 0.2 μm filter before lyophilizing

This material was prepared from bovine plasma obtained from cattle of less than 30 months of age and that have received ante- and post-mortem health inspections from APHIS Inspectors within facilities inspected by the USDA.



Quality Assurance