

Product: **Bovine Serum Albumin (BSA)**
Fraction V—Low Electrolyte
 Powder

Catalog #: 700-109P

Lot #: C13H09H

Origin: United States (100%)

Manufacture Date: Feb 2016

Expiration Date: Mar 2019

<i>TEST</i>	<i>METHODOLOGY</i>	<i>SPECIFICATION</i>	<i>RESULTS</i>
Microbiological Testing			
Microbial Growth		Not Detected	Not Detected
Viral Testing	Modified 9 CFR 113.53c		
Extraneous Viruses (FAb)			
Bluetongue Virus		Not Detected	Not Detected
Vesicular Stomatitis Virus		Not Detected	Not Detected
Physical Testing			
pH @ RT, 10% Solution		4.9-5.2	5.1
Conductivity, 10% Solution		≤2,200 μS	530 μS
Biochemical Testing			
Purity	HPLC	≥98%	100%

Filtered through a 0.2 μm filter before lyophilizing

This material was prepared from bovine plasma obtained from cattle of less than 30 months of age and that have received ante- and post-mortem health inspections from APHIS Inspectors within facilities inspected by the USDA.



 Created by



 Verified by