

Product: **Bovine Serum Albumin (BSA)**
Fraction V—RIA Grade
 Powder


Catalog #: 700-103P
Lot #: C32H09H
Origin: United States (100%)

Manufacture Date: April 2016
Expiration Date: April 2019


TEST	METHODOLOGY	SPECIFICATION	RESULTS
Microbiological Testing			
Microbial Growth		≤3 cfu/mL	No Growth
Viral Testing			
Extraneous Viruses (FAb)			
Bluetongue Virus		Negative	Negative
Vesicular Stomatitis Virus		Negative	Negative
Physical Testing			
pH @ RT, 10% Solution		6.7-7.3	6.8
Conductivity, 10% Solution		≤3,000 μS	1100 μS
Biochemical Testing			
Purity	HPLC	≥98%	100%
Protease		Not detected	Not detected
IgG		Below detectable limits	Below detectable limits

Filtered through a 0.2 μm filter before lyophilizing

This material was prepared from bovine plasma obtained from cattle of less than 30 months of age and that have received ante- and post-mortem health inspections from APHIS Inspectors within facilities inspected by the USDA.



 Created by



 Verified by