



GEMINI
BIO-PRODUCTS

Certificate of Analysis

Created April 10, 2018

Product: **Bovine Serum Albumin (BSA)**
Fraction V—Cohn
Powder

Catalog #: 700-108P

Lot #: C67H17J

Origin: United States (100%)

Manufacture Date: Mar 2018

Expiration Date: Mar 2023

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Microbiological Testing			
Microbial Growth	Pour Plate, 30-35°C, 48 hrs	Passed	Passed
Viral Testing			
Bluetongue Virus	9 CFR 113	Not Detected	Not Detected
Vesicular Stomatitis Virus	VN VS NJ/IND	Not Detected	Not Detected
Physical Testing			
pH	pH Meter	Report	7.00
Biochemical Testing			
Purity	Protein Electrophoresis	≥98%	98.2%

This material was prepared from bovine plasma obtained from cattle that have received ante- and post-mortem health inspections from APHIS Inspectors within facilities inspected by the USDA.


Created by


Verified by