

Product: **Bovine Serum Albumin (BSA)**
Fraction V - Standard Grade
 Powder

Catalog #: 700-100P
Lot #: C77G09H
Origin: United States (100%)

Manufacture Date: July 2016
Expiration Date: July 2019

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Viral Testing			
Extraneous Viruses (FAb)			
Bluetongue Virus		Negative	Negative
Vesicular Stomatitis Virus		Negative	Negative
Physical Testing			
pH @ RT, 10% Solution		6.7-7.3	6.8
Conductivity, 10% Solution		≤3,000 μS	1000 μS
Biochemical Testing			
Purity	HPLC	≥98%	100%

This material was prepared from bovine plasma obtained from cattle that have received ante- and post-mortem health inspections from APHIS Inspectors within facilities inspected by the USDA.



 Quality Assurance