



Certificate of Analysis

Created October 6, 2016

Product: **Bovine Serum Albumin (BSA)**
Fraction V - Standard Grade
Powder

Catalog #: 700-100P
Lot #: C83G09H
Origin: United States (100%)

Manufacture Date: Aug 2016
Expiration Date: Aug 2019

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Viral Testing			
Extraneous Viruses (FAb)			
Bluetongue Virus		Negative	Negative
Vesicular Stomatitis Virus		Negative	Negative
Physical Testing			
pH @ RT, 10% Solution		6.7-7.3	7.0
Conductivity, 10% Solution		≤3,000 μS	1300 μS
Biochemical Testing			
Purity	HPLC	≥98%	100%

This material was prepared from bovine plasma obtained from cattle that have received ante- and post-mortem health inspections from APHIS Inspectors within facilities inspected by the USDA.


Quality Assurance