



Certificate of Analysis

Created February 26, 2016

Product: **Donor Horse Serum**
Dialyzed
0.2µm sterile-filtered

Catalog #: 100-158

Lot #: E13H03H

Origin: United States (100%)

Manufacture Date: Feb 2016

Expiration Date: Feb 2021

DATA SHEET

Description

Material derived from donor horse serum, dialyzed at 10,000 mw cut-off, then 0.2µm filtered.

All material should be handled as though capable of transmitting infectious agents.



QUALITY CONTROL