



**GEMINI**  
BIO-PRODUCTS

# Certificate of Analysis

Created February 12, 2021

**Product:** **Trypsin: EDTA Solution**  
0.25% trypsin—0.53mM EDTA  
Sterile-filtered

**Catalog #: 400-151**  
**Lot #: F021008**

**Manufacture Date: Jan 2021**  
**Expiration Date: Jan 2022**

| TEST                           | METHODOLOGY | SPECIFICATION | RESULTS     |
|--------------------------------|-------------|---------------|-------------|
| <b>Microbiological Testing</b> |             |               |             |
| Sterility                      | USP<71>     |               |             |
| Bacteria                       |             | No Growth     | No Growth   |
| Fungi                          |             | No Growth     | No Growth   |
| <b>Physical Testing</b>        |             |               |             |
| Osmolality                     | USP<785>    | Report        | 282 mOsm/kg |
| pH @ RT                        | USP<791>    | Report        | 7.74        |

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.*

   *Amu* 12 FEB 2021

Created by/ date

   *Chelin B. Hut* 12 FEB 2021

Verified by/ date