



**Product:** **Polymixin B Solution (50mg/mL)**  
Sterile-filtered

**Catalog #: 400-177**  
**Lot #: F68T00L**

**Manufacture Date: May 2020**  
**Expiration Date: May 2022**

TEST	METHODOLOGY	SPECIFICATION	RESULTS
<b>Microbiological Testing</b>			
Sterility	USP<71>		
Bacteria		Not Detected	Not Detected
Fungi		Not Detected	Not Detected
<b>Physical Testing</b>			
Osmolality	USP<785>	Report	53 mOsm/kg
pH @ RT	USP<791>	Report	5.81
<b>Active Unit Activity</b>			
Polymixin B Sulfate Powder			9680 IU/mg

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.*

*[Handwritten Signature]* 18 May 2020

Created by/ date

*[Handwritten Signature]* 18 May 2020

Verified by/ date