



Product: **L-Glutamine (Stabilized)**
C₈H₁₅N₃O₄
Powder

Catalog #: 400-106P
Lot #: F75R14I

Manufacture Date: Jul 2017
Expiration Date: Jul 2019

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Biological Testing Assay (Dried Basis)		≥ 98.5%	99.6%
Biochemical Testing Chloride		≤ 0.02%	< 0.02%
Heavy Metals		≤ 10 ppm	< 10 ppm



Created by



Verified by