



Product: **L-Glutamine Solution (Stabilized)**
L-alanyl-L-glutamine
Sterile-filtered

Catalog #: 400-106
Lot #: F75T00L

Manufacture Date: Aug 2020
Expiration Date: Aug 2022

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Microbiological Testing			
Sterility	USP<71>		
Bacteria		Not Detected	Not Detected
Fungi		Not Detected	Not Detected
Biochemical Testing			
Osmolality	USP<785>	Report	464 mOsm/kg
pH	USP<791>	Report	5.83

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.

[Signature] 21 Aug 2020

Created by/ date

[Signature] 21 AUG 2020

Verified by/ date