



GEMINI
BIO-PRODUCTS

Certificate of Analysis

Created April 3, 2019

Product: mLIF (Leukemia Inhibitory Factor-mouse)

Catalog #: 400-495
Lot #: F90S00K

Manufacture Date: Feb 2019
Expiration Date: Feb 2021

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Biological Testing			
Endotoxin		Report	4.41 EU/ug
Microbiological Testing			
Sterility			
Bacteria		Not Detected	Not Detected
Fungi		Not Detected	Not Detected

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.

Emelyne D. Davis 03 Apr 2019
Created by/ date

APM MA 03 Apr 2019
Verified by/ date