



**Product:** **Recombinant Human/ Mouse FGF-8**  
**(Fibroblast growth factor 8)**  
Powder

**Catalog #: 300-219P**  
**Lot #: G32L40I**  
**Assayed: Nov 2017**

**Date Received:** \_\_\_\_\_

**Expiration Date:\*** \_\_\_\_\_

TEST	METHODOLOGY	RESULTS
Biological Testing		
Endotoxin	LAL	<0.05 EU/μg
Biochemical Testing		
Purity	SDS PAGE	>95%
Biological Activity	NR6-R 3T3 cell proliferation	1.1 x 10 <sup>6</sup> U/mg

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.*

*Erlyne Barber* 02 Apr 2019  
Created by/ date

*[Signature]* 02 Apr 2019  
Verified by/ date

\*Expiration Date is one year from the date of receipt when stored at -20 to -80°C.