



**Product:** **Recombinant Human IL-17AF**  
Interleukin 17AF  
Powder

**Catalog #: 300-141P**  
**Lot #: G63L40Z**  
**Assayed: May 2009**

**Date Received:** \_\_\_\_\_  
**Expiration Date:\*** \_\_\_\_\_

TEST	METHODOLOGY	RESULTS
Biological Testing		
Endotoxin	LAL	0.31 EUs/ $\mu$ g
Biochemical Testing		
Purity	SDS PAGE	>95%
Biological Activity ED <sub>50</sub>	Production of IL-6 from mouse 3T3 cells	3.1 x 10 <sup>5</sup> U/mg

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.*

Oliver 23 sep 2019  
Created by/ date

MJM 23 sep 2019  
Verified by/ date

\*Expiration Date is one year from the date of receipt when stored at -20 to -80°C.