



Product: **Human TGF-beta-3**
Human Recombinant
Solution

Catalog #: 300-192
Lot #: G67K40E

Assayed: Feb 2013

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Biological Testing Endotoxin	LAL	≤1 EU/μg	0.078 EU/μg
Biochemical Testing Purity	SDS PAGE/RP-HPLC	≥95%	>95%
Biological Activity ED ₅₀	Inhibition of IL-4-induced HT-2 cell proliferation	≥1.0x10 ⁶ U/mg	2.6x10 ⁷ U/mg

Created by

Verified by

Expiration Date

One year from the date of receipt when stored at 4°C as supplied.

Date Received: _____

Expiration Date: _____