



GEMINI
BIO-PRODUCTS

Certificate of Analysis

Created February 22, 2018

Product: **Recombinant Human IL-35 Heterodimer**
(Interleukin-35)
Powder

Catalog #: 300-212P
Lot #: G68K40E
Assayed: Jan 2013

Date Received: _____

Expiration Date:* _____

TEST	METHODOLOGY	RESULTS
Biological Testing Endotoxin	LAL	1.69 EU/ μ g
Biochemical Testing Purity	SDS PAGE	>95%

Created by

Verified by

*Expiration Date is one year from the date of receipt when stored at -20°C to -80°C as supplied.