



Product: **Recombinant Mouse IL-11**
(Interleukin-11)
Powder

Catalog #: 300-313P

Lot #: G75K40I

Assayed: Dec 2017

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Biological Testing Endotoxin	LAL	≤1 EU/μg	<0.052 EU/μg
Biochemical Testing Purity	SDS-PAGE	≥95%	>95%
Biological Activity ED ₅₀	B9 Cell Proliferation	≥4.0 x 10 ³ U/mg	7.4 x10 ³ U/mg



Created by



Verified by

Expiration Date

One year from the date of receipt when stored at -20 to -80 °C as supplied.

Date Received: _____

Expiration Date: _____