

**Product:** **Human Serum Albumin**  
**Fatty Acid Free**  
 Powder

**Catalog #: 800-107P**

**Lot #: H11R09G**

**Origin: United States (100%)**

**Manufacture Date: Oct 2015**

**Expiration Date: Sep 2018**

TEST	METHODOLOGY	SPECIFICATION	RESULTS
<b>Microbiological Testing</b>			
Microbial Growth		<2.99 cfu/mL	0 cfu/mL
<b>Viral Testing</b>			
HBsAg		Non-Reactive	Non-Reactive
Anti-HCV		Non-Reactive	Non-Reactive
Anti-HIV-1/HIV-2		Non-Reactive	Non-Reactive
HIV-1/HIV-2-Ab		Not Detected	Not Detected
HCV-Ab		Not Detected	Not Detected
<b>Physical Testing</b>			
pH @ RT			6.9
<b>Biochemical Testing</b>			
Fatty Acid			<0.01%
Purity	HPLC		96%

*Results shown were obtained by carefully performed methods believed to be reliable. However, since some results may vary for specific tests depending upon methodology and other variables, it is suggested that tests for which results are particularly important be repeated by the user of this product.*

*All human plasma is collected from stringently screened donors at FDA-licensed collection centers in the United States. All donors are tested and found to be negative for Indirect Antiglobulin.*



Quality Assurance