



# Certificate of Analysis

Created October 24, 2016

**Product:** **Human Serum AB (Off the Clot)**  
**Heat Inactivated**  
 0.1µm sterile-filtered

**Catalog #: 100-318**  
**Lot #: H47R00H**  
**Origin: United States (100%)**

**Manufacture Date: Aug 2016**  
**Expiration Date: Aug 2021**

TEST	METHODOLOGY	SPECIFICATION	RESULTS
<b>Biological Testing</b>			
Endotoxin	USP<85>		2.48 EU/mL
Hemoglobin		Report	2.55 mg/dL
<b>Microbiological Testing</b>			
Sterility			
Bacteria		Not Detected	Not Detected
Fungi		Not Detected	Not Detected
Mycoplasma		Not Detected	Not Detected
<b>Viral Testing Donor</b>			
HBsAg		Non-Reactive	Non-Reactive
HBV		Non-Reactive	Non-Reactive
HCV		Non-Reactive	Non-Reactive
HIV-1/HIV-2		Non-Reactive	Non-Reactive
Syphilis		Negative	Negative
HIV-1/HIV-2-NAT		Not Detected	Not Detected
HCV-NAT		Not Detected	Not Detected
<b>Physical Testing</b>			
Osmolality		260-350 mOsm/kg	291 mOsm/kg
pH @ RT		Report	7.63

**Product:** **Human Serum AB (Off the Clot)**  
**Heat Inactivated**  
 0.1µm sterile-filtered

**Catalog #: 100-318**

**Lot #: H47R00H**

**Origin: United States (100%)**

**Manufacture Date: Aug 2016**

**Expiration Date: Aug 2021**

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Albumin			4.2 g/dL
ALT (SGPT)			15 U/L
AST (SGOT)			15 U/L
Bilirubin, Total			0.2 mg/dL
BUN			14 mg/dL
Calcium			2.33 mM
Chloride			104 mM
Cholesterol			158 mg/dL
Creatinine			1.04 mg/dL
Glucose			5.61 mM
Phosphorus			1.20 mM
Potassium			4.5 mM
Protein, Total			7.2 g/dL
Sodium			141 mM
Triglycerides			106 mg/dL
Uric Acid			5.6 mg/dL

*Results shown were obtained by carefully performed methods believed to be reliable. However, since some results may vary for specific tests depending upon methodology and other variables, it is suggested that tests for which results are particularly important be repeated by the user of this product.*

*All human blood products are collected from stringently screened donors at FDA-licensed collection centers in the United States. All donors are tested and found to be negative for Indirect Antiglobulin.*

*Precipitates may develop in this product upon freezing and/or thawing; this occurrence does not impact culture performance.*

*For further manufacturing or laboratory use. Not for human or animal consumption.*



Quality Assurance