



Certificate of Analysis

Created December 7, 2020

Product: **Normal Human Serum**
Individual Male Donor
Off the Clot

Catalog #: 900-110M
Lot #: H71Y03L

Draw Date: Sep 2020


DATA SHEET

Gender Male

Testing This material has been tested by an FDA approved test and found to be negative for Chagas, HBsAg, HIV-1/2, HIV-1 NAT, HBV-NAT, HCV-NAT, HCV, WNV, and Syphilis.

All human materials should be handled as though capable of transmitting infectious agents. For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.

 07DEC2020
Created by/ date

 07DEC-2020
Verified by/ date