



Certificate of Analysis

Created November 9, 2016

Product: **Grace's Insect Medium**
With L-Glutamine
Sterile-filtered

Catalog #: 600-310
Lot #: M54Q20H

Manufacture Date: May 2016
Expiration Date: Nov 2017

<i>TEST</i>	<i>METHODOLOGY</i>	<i>SPECIFICATION</i>	<i>RESULTS</i>
Biological Testing			
Endotoxin	USP<85>	Report	<0.25 EU/mL
Microbiological Testing			
Sterility			
Bacteria		Not Detected	Not Detected
Fungi		Not Detected	Not Detected
Physical Testing			
Osmolality			362 mOsm/kg
pH @ RT			6.30



Quality Assurance