

Product: **TNM-FH Insect Medium**
 (Supplemented Grace's)
 With 3330 mg/L Yeastolate and 3330 mg/L Lactalbumin Hydrolysate
 Sterile-filtered

Catalog #: 600-311
Lot #: M61Q20H

Manufacture Date: Dec 2016
Expiration Date: Jun 2018

<i>TEST</i>	<i>METHODOLOGY</i>	<i>SPECIFICATION</i>	<i>RESULTS</i>
Biological Testing			
Endotoxin	USP<85>	Report	<1 EU/mL
Microbiological Testing			
Sterility			
Bacteria		Not Detected	Not Detected
Fungi		Not Detected	Not Detected
Physical Testing			
Osmolality			382 mOsm/kg
pH @ RT		6.0 - 6.4	6.03



 Quality Assurance