

Product: **COHSII Media**
 0.2µm sterile-filtered

Catalog #: 900-750
Lot #: M68Q00I

Manufacture Date: May 2017
Expiration Date: May 2019

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Biological Testing			
Endotoxin	USP <85>	Report	<0.5 EU/mL
Microbiological Testing			
Sterility			
Bacteria		Not Detected	Not Detected
Fungi		Not Detected	Not Detected
Mycoplasma	Barile, MF and Kern, J (1971)	Not Detected	Not Detected
General Safety Testing*	USP <88>	Pass	Pass
Physical Testing			
Osmolality		275 - 340	304 mOsm/kg
pH		7.3 - 7.5	7.36

*Note: The sample meets requirements of the General Safety Test under the United States Pharmacopeia and National Formulary, Chapter <88> Biological Reactivity Tests.



 Created by



 Verified by