




Product: **Hank's Balanced Salt Solution**
Sterile-filtered

Catalog #: 600-101
Lot #: M58S00L


Manufacture Date: Apr 2020
Expiration Date: Apr 2022

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Biological Testing			
Endotoxin	USP <85>	Report	<0.500 EU/mL
Microbiological Testing			
Sterility	USP<71>		
Bacteria		Not Detected	Not Detected
Fungi		Not Detected	Not Detected
Physical Testing			
Osmolality		270-305 mOsm/Kg	279 mOsm/kg
pH		6.7-7.8	7.8

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.



 Created by/ date



 Verified by/ date